


X-tension


Instruksjoner for bruk

Dette produktet er designet til å motstå eller stoppe dorsifleksjon med Ankel-Fot-Ortoser (AFO) .
X-tension is designed to resist or stop dorsiflexion.



- 1) Hvis din ortose krever et stopp, tilpasses din AFO sagittalt før du installerer X-tension. (SNAPstop er det eneste plantarfleksjonsstopp, som passer til X-tensions geometri)
If your orthosis requires a stop, sagittally align your AFO before installing X-tension. (SNAPstop is the only plantar flexion stop that fits the X-tension geometry)
- 2) Plasser X-tension-båndet på den bakerste delen av AFOen med SNAPstop-vinduet direkte over SNAPstop. (Hvis SNAPstop ikke anvendes, plasseres bunnen av SNAPstop-vinduet på nivå med akselen i ankelleddet)
Position the X-tension band on the posterior aspect of the AFO with the SNAPstop Window directly over the SNAPstop. (If SNAPstop is not used, position the bottom of the SNAPstop Window level with the axis of the ankle joint).
- 3) Markér plasseringen av de to nederste monteringshullene på plasten med en markeringsblyant eller brenner.
Mark the location of the two lower mounting holes on the plastic with a marking pencil or with hot poker.
- 4) Bor et 1/8" hul (3 mm) på de markerte stedene.
Drill a 3 mm hole (1/8") hole at the marked points.
- 5) Fastgjør det nedre båndet med røniagler (for dynamisk X-tension) eller kobbermagler (for statisk X-tension)
Secure the lower band with Speedy rivets (for Dynamic X-tension) or copper rivets (for Static X-tension).
- 6) Når de to bunnhullene er fastgjort, strekkes X-tension-båndet over AFO en og de to øverste monteringshullene markeres.
Once the two bottom holes are secured, stretch the X-tension band over the AFO and mark the two upper mounting holes.
- 7) Bor et 1/8" (3 mm) hul på de markerte stedene. For at øke spenningen i X-tension-båndet, bores hullene litt oppover og utover fra merkene.
Drill a 3 mm (1/8") hole at the marked points. To increase tension in the X-tension band, drill your holes slightly up and out from your marks.

- 8) Strekkes X-tension-båndet (hvis dynamisk) og juster hullene. Fastgjør det øverste båndet med rømiagler (for dynamisk X-tension) eller kobbernagler (for statisk X-tension).
Stretch the X-tension band slightly (if Dynamic) and align the holes. Secure the upper band with Speedy rivets (for Dynamic X-tension) or copper rivets (for Static X-tension).



Dette produktet bør kun anvendes til dets opprinnelige formål og bør fremstilles under veiledning av en sertifisert bandasjist i henhold til de medfølgende instruksjonene

This product should only be used for its intended purpose and should be fabricated under the guidance of a certified Orthotist per the provided instructions.



www.launchpad-europe.eu

Manufactured by LaunchPad O&P
2545 Chicago Ave, Suite 507, Minneapolis, MN 55404 USA

Authorized Representative EU:
Dirame Ortho, ASSESTEENWEG 27-29, B-1740 TERNAT, Belgium
T: +32 (0)2 582 82 50 / F: +32 (0)2 582 61 12

